

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA808825

RPS6KC1 Mouse Monoclonal Antibody [Clone ID: OTI10F11]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI10F11
Applications: IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 1-305 of human

RPS6KC1(NP_036556) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 118.5 kDa

Gene Name: ribosomal protein S6 kinase C1

Database Link: NP 036556 Entrez Gene 26750 Human

Background: May be involved in transmitting sphingosine-1 phosphate (SPP)-mediated signaling into the


cell. [UniProtKB/Swiss-Prot Function]

Synonyms: humS6PKh1; RPK118; S6K-delta-1; S6PKh1


Protein Families: Druggable Genome, Protein Kinase



Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY RPS6KC1 ([RC211142], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RPS6KC1 (1:2000). Positive lysates [LY402211] (100ug) and [LC402211] (20ug) can be purchased separately from OriGene.



Immunohistochemical staining of paraffinembedded Human Kidney tissue using anti-RPS6KC1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA808825) (1:150)