

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA809094

GATA3 Mouse Monoclonal Antibody [Clone ID: OTI5C11]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI5C11

Applications: FC, IF, IHC, WB

Recommend Dilution: IHC 1:200

Reactivity: Human

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 155-443 of human

GATA3 (NP 002042) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 47.9 kDa

Gene Name: GATA binding protein 3

Database Link: NP 001002295 Entrez Gene 2625 Human

Background: This gene encodes a protein which belongs to the GATA family of transcription factors. The

protein contains two GATA-type zinc fingers and is an important regulator of T-cell

development and plays an important role in endothelial cell biology. Defects in this gene are the cause of hypoparathyroidism with sensorineural deafness and renal dysplasia. [provided

by RefSeq, Nov 2009]

Synonyms: HDR; HDRS


Protein Families: Adult stem cells, ES Cell Differentiation/IPS, Stem cell relevant signaling - JAK/STAT signaling

pathway, Transcription Factors






Product images:




Western blot analysis of extracts (35ug) from 3 different cell lines by using anti-GATA3 monoclonal antibody (1:200).




Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-GATA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA809094) (1:200)






Immunohistochemical staining of paraffinembedded Carcinoma of Human kidney tissue using anti-GATA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA809094) (1:200)



Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-GATA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA809094) (1:200)




Immunohistochemical staining of paraffinembedded Carcinoma of Human bladder tissue using anti-GATA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA809094) (1:200)




MCF7


HUVEC





Immunocytochemistry staining of MCF-7 cells using anti-GATA3 mouse monoclonal antibody (TA809094) (Left). The right is negative control. (1:2000)





Flow cytometric Analysis of MCF-7 cells, using anti-GATA3 antibody (TA809094), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).The right is HUVEC cells as negtive control (1:100).