

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA809101**

## **ZNF230 Mouse Monoclonal Antibody [Clone ID: OTI8F8]**

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI8F8
Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human

Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Full length human recombinant protein of human ZNF230 (NP\_006291) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 54.4 kDa

**Gene Name:** zinc finger protein 230

Database Link: NP 006291 Entrez Gene 7773 Human

**Background:** May be involved in transcriptional regulation. [UniProtKB/Swiss-Prot Function]


Synonyms: FDZF2

**Protein Families:** Transcription Factors





## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ZNF230 ([RC208201], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZNF230 (1:2000). Positive lysates [LY416742] (100ug) and [LC416742] (20ug) can be purchased separately from OriGene.