

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA809330

C1orf83 (TCEANC2) Mouse Monoclonal Antibody [Clone ID: OTI9A4]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI9A4

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human
Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human TCEANC2(NP_694580) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 24 kDa

Gene Name: transcription elongation factor A N-terminal and central domain containing 2

Database Link: NP 694580 Entrez Gene 127428 Human


Synonyms: C1orf83

Protein Families: Transcription Factors





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TCEANC2 ([RC206749], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TCEANC2 (1:2000). Positive lysates [LY407172] (100ug) and [LC407172] (20ug) can be purchased separately from OriGene.