

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA809736**

### CD161 (KLRB1) Mouse Monoclonal Antibody [Clone ID: OTI1D8]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI1D8
Applications: IHC, WB

**Recommend Dilution:** WB 1:500, IHC 1:500

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 67-225 of human

KLRB1(NP\_002249) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 25.2 kDa

**Gene Name:** killer cell lectin like receptor B1

Database Link: NP 002249 Entrez Gene 3820 Human

**Background:** Natural killer (NK) cells are lymphocytes that mediate cytotoxicity and secrete cytokines after

immune stimulation. Several genes of the C-type lectin superfamily, including the rodent NKRP1 family of glycoproteins, are expressed by NK cells and may be involved in the regulation of NK cell function. The KLRB1 protein contains an extracellular domain with several motifs characteristic of C-type lectins, a transmembrane domain, and a cytoplasmic domain. The KLRB1 protein is classified as a type II membrane protein because it has an

external C terminus. [provided by RefSeq, Jul 2008]


Synonyms: CD161; CLEC5B; hNKR-P1A; NKR; NKR-P1; NKR-P1A; NKRP1A

**Protein Families:** Transmembrane






## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY KLRB1 ([RC216459], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-KLRB1 (1:500). Positive lysates [LY419438] (100ug) and [LC419438] (20ug) can be purchased separately from OriGene.



Immunohistochemical staining of paraffinembedded Human spleen tissue within the normal limits using anti-KLRB1 mouse monoclonal antibody. (Heat-induced epitope retrieval by Tris-EDTA (1:500)