

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA809737**

### TRIP13 Mouse Monoclonal Antibody [Clone ID: OTI2F5]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI2F5
Applications: IHC, WB

Recommend Dilution: WB 1:500~2000, IHC 1:2000

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 31-148 of human

TRIP13(NP\_001159732) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Gene Name: thyroid hormone receptor interactor 13

Database Link: NP 001159732 Entrez Gene 9319 Human

**Background:** This gene encodes a protein that interacts with thyroid hormone receptors, also known as

hormone-dependent transcription factors. The gene product interacts specifically with the ligand binding domain. This gene is one of several that may play a role in early-stage non-

small cell lung cancer. [provided by RefSeq, Oct 2009]


Synonyms: 16E1BP

**Protein Families:** Druggable Genome, Stem cell - Pluripotency, Transcription Factors






## **Product images:**




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TRIP13 ([RC228224], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TRIP13 (1:2000).




Western blot analysis of extracts (35ug) from 3 different cell lines by using anti-TRIP13 monoclonal antibody (293T: human; HeLa: human; Jurkat: human; ) (1:500).






Equivalent amounts of cell lysates (10 ug per lane) of wild-type 293T cells (WT, Cat# LC810293T) and TRIP13-Knockout 293T cells (KO, Cat# [LC811410]) were separated by SDS-PAGE and immunoblotted with anti-TRIP13 monoclonal antibody TA809737, (1:500). Then the blotted membrane was stripped and reprobed with anti-HSP90AB1 antibody ([TA500494]) as a loading control.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-TRIP13 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA809737) (1:2000)



Immunohistochemical staining of paraffinembedded Human endometrium tissue within the normal limits using anti-TRIP13 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA809737) (1:2000)





Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-TRIP13 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA809737) (1:2000)