

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA809840

CBR1 Mouse Monoclonal Antibody [Clone ID: OTI4B4]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4B4
Applications: IHC, WB

Recommend Dilution: WB 1:500~2000, IHC 1:2000

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human CBR1 (NP_001748) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 30.2 kDa

Gene Name: carbonyl reductase 1

Database Link: NP 001748 Entrez Gene 873 Human

Background: The protein encoded by this gene belongs to the short-chain dehydrogenases/reductases

(SDR) family, which function as NADPH-dependent oxidoreductases having wide specificity for

carbonyl compounds, such as quinones, prostaglandins, and various xenobiotics.

Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq,

Nov 2013]

Synonyms: CBR; hCBR1; SDR21C1


Protein Families: Druggable Genome

Protein Pathways: Arachidonic acid metabolism, Metabolic pathways






Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CBR1 ([RC204950], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CBR1 (1:2000). Positive lysates [LY400669] (100ug) and [LC400669] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 4 different tissue lysates by using anti-CBR1 monoclonal antibody (1:500).






Western blot analysis of extracts (35ug) from 4 different cell lines by using anti-CBR1 monoclonal antibody (1:500).



Equivalent amounts of cell lysates (10 ug per lane) of wild-type A549 cells and CBR1-Knockout A549 cells (KO, Cat# [LC800036]) were separated by SDS-PAGE and immunoblotted with anti-CBR1 monoclonal antibody TA809840 (1:500). Then the blotted membrane was stripped and reprobed with anti-HSP90 antibody as a loading control.



Immunohistochemical staining of paraffinembedded Carcinoma of Human liver tissue using anti-CBR1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA809840) (1:2000)