

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA809996

RHOH Mouse Monoclonal Antibody [Clone ID: OTI9F4]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI9F4
Applications: WB

December of Bile Control NAD 4-200

Reactivity: WB 1:2000
Host: WB 1:2000
WB 1:2000
Human
Mouse

Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 97-190 of human

RHOH(NP_004301) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 21.1 kDa

Gene Name: ras homolog family member H

Database Link: NP 004301 Entrez Gene 399 Human

Background: The protein encoded by this gene is a member of the Ras superfamily of guanosine

triphosphate (GTP)-metabolizing enzymes. The encoded protein is expressed in

hematopoietic cells, where it functions as a negative regulator of cell growth and survival.

This gene may be hypermutated or misexpressed in leukemias and lymphomas.

Chromosomal translocations in non-Hodgkin's lymphoma occur between this locus and B-cell CLL/lymphoma 6 (BCL6) on chromosome 3, leading to the production of fusion transcripts. Alternative splicing in the 5' untranslated region results in multiple transcript variants that

encode the same protein. [provided by RefSeq, May 2013]

Synonyms: ARHH; TTF


Protein Families: Transcription Factors

Protein Pathways: Leukocyte transendothelial migration





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY RHOH ([RC204324], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RHOH (1:2000). Positive lysates [LY401372] (100ug) and [LC401372] (20ug) can be purchased separately from OriGene.